HIGH PERFORMANCE RESEARCH COMPUTING ACES: Introduction to CryoSPARC for Cryo-EM Data Processing in Collaboration with the Laboratory for Biomolecular Structure and Dynamics



High Performance Research Computing DIVISION OF RESEARCH

Fall 2023

CryoSPARC on ACES


- CryoSPARC Academic License ID
- Resources and Limitations
- CryoSPARC on the ACES Portal
- Group data directories
- CryoSPARC Tutorial

CryoSPARC Academic License ID

- CryoSPARC provides a personal Academic License ID free of charge for TAMU staff and students
- A CryoSPARC Academic License ID is required to launch the CryoSPARC ACES portal app
- Use your academic email address for the Academic License
 - <u>https://cryosparc.com/download</u>

Resources for Running CryoSPARC


- CryoSPARC is available as an Interactive App on the ACES Portal
 - <u>https://portal-aces.hprc.tamu.edu</u>
- ACES hardware overview (can change since ACES is a composable resource cluster)
 - 30 x H100 GPUs (2 x H100 GPUs per node)
 - 6 x A30 GPUs (6 x A30 GPUs per node)
 - 70 x CPU only nodes
 - 10 x other specialized nodes not for CryoSPARC
- Specify enough time to allow your processing to complete
- If you launch a job for 24 hours and you finish your work in 12 hours, delete your job to make GPUs available for other jobs



Resource Limitations

- Can only launch one CryoSPARC portal session at a time
 can run multiple CryoSPARC jobs within a portal session
- CPUs
 - Some CryoSPARC tasks do not require a GPU such as importing images
- GPUs
 - The GPU queue can get busy at times and it may take over an hour for your job to launch
 - Cancel your pending job if you will be away from your computer for a long time and it hasn't started yet
- Submitted jobs can have delayed start times due to unavailability of GPUs or scheduled system maintenance

Accessing the ACES Portal




HPRC webpage: hprc.tamu.edu

ĂМ


Accessing ACES via the ACCESS Portal

Log-in using your ACCESS credentials.



ACES CryoSPARC Portal App

portal-aces.hprc.tamu.edu



CryoSPARC

This app will launch CryoSPARC on one ACES cluster compute node.

Cryo	SPARC license id (required)
	You must get an individual license id from
	cryosparc.com/download
	After you click the blue Launch button below, it may
	take up to 5 minutes before CryoSPARC is available
	Use the following to login to your CryoSPARC session
	 username: admin@cryo.edu
	o password: admin
	You can only launch one CryoSPARC portal app job
	but you can run multiple computational jobs within
	CryoSPARC
Num	ber of hours (max 168)
3	\$

8

Launch CryoSPARC the first time using 'CPU only' in order to initialize the database.


Shell Access via the ACES Portal



Display GPU Availability on ACES

See GPU configuration and current availability using the ACES portal shell access

- <u>https://portal-aces.hprc.tamu.edu</u>
- select "Clusters" then ">_aces Shell Access"
- at the command line prompt, enter: gpuavail



AVAILABILITY AVAILABILITY						
NODE NAME	GPU TYPE	GPU COUNT	GPU AVAIL	CPU AVAIL	GB MEM AVAIL	
ac027	h100	4	1	72	248	
ac048	h100	2	2	96	488	
ac049	h100	2	1	88	408	
ac050	h100	2	2	96	488	
ac064	a30	6	6	96	488	

ACES Maintenance

- ACES will be offline during routinely scheduled maintenances which may occur a few times each year
- Submitted jobs will remain PENDING and not run if the job script wall time overlaps with the maintenance start time
- When maintenance is complete, scheduled jobs will automatically begin running based on priority
- Use the **maintenance** command to see if a maintenance is scheduled and to display the time until the start of the maintenance


maintenance

The scheduled 35 hour ACES maintenance will start in:

3 days 21 hours 32 minutes


Scheduled jobs will not start if they overlap with this maintenance window.

ACES CryoSPARC Running Job



ACES CryoSPARC Portal App

portal-aces.hprc.tamu.edu



Delete your CryoSPARC "CPU only" job and launch a new CryoSPARC job using 1 x H100 GPU for the tutorial


CryoSPARC This app will launch CryoSPARC on one ACES cluster compute node. Version 4.2.1 CryoSPARC license id (required) You must get an individual license id from cryosparc.com/download After you click the blue Launch button below, it may take up to 5 minutes before CryoSPARC is available Use the following to login to your CryoSPARC session: username: admin@cryo.edu password: admin · You can only launch one CryoSPARC portal app job but you can run multiple computational jobs within CryoSPARC Number of hours (max 96) \$ 3 Node type H100 GPU Number of GPUs \$

CryoSPARC Visualization of Resource Utilization

Review CPU and GPU usage for a Job

The **jobstats** utility automatically runs for each CryoSPARC portal job and monitors CPU and GPU resource usage and creates a graph for each.

- CPU stats monitors all cores regardless of how many were configured for the job
 - CryoSPARC jobs use a fraction of total cores and memory if selecting one of many GPUs on a compute node to allow others to use the additional GPUs on the compute node
- GPU stats will create a graph for each GPU that was configured for the job




Review CPU and GPU usage for a CryoSPARC Job

- At the bottom of the CryoSPARC portal app, access the CryoSPARC data root directory
- 2. Click 'output'

ĀΜ

- 3. Click 'Modified at' to sort newest session on top
- 4. Click the Name of the session session you want to view based on the date
- 5. Click the .png image files to view the CPU and GPU usage stats



Group Data Directories

- Group data directories can be used to share input files or to have a group work space
- Send a request to the HPRC helpdesk to create a group
 - <u>help@hprc.tamu.edu</u>
 - provide a group name and usernames of members
 - group directories have their own disk quotas separate from individual users
- You will need to mount your group directory when launching the CryoSPARC portal app

Optional group directory to mount

- You must already be a member of the group
- Example values:
 - /junjiez
 - /scratch/group/davislab

CryoSPARC Tutorial